

AUTOMATED MRM/PRM QUANTITATION · VALIDATION BENCHMARK

Fully automated analysis of MRM/PRM data

Targeted MRM/PRM is the workhorse of biomarker validation and clinical assays. The bottleneck isn't the mass spec — it's turning thousands of chromatograms into trustworthy numbers, one manually-reviewed peak at a time.

The problem

Manual peak review is slow, subjective, and doesn't scale. A single validation run can tie up an experienced analyst for hours to days — and the answer still depends on who did the picking.

What MRMpipe does

From raw vendor files to a validated quant table, fully automated: peak detection, integration, calibration and QC — with light/heavy co-elution enforced at pick time, so the machine doesn't make the mistakes a rushed reviewer would.

The only question that matters

Does automation actually match an expert's curated result? We put it to the test on a real **anti-doping panel** — **918 samples, 35 peptides** — benchmarking MRMpipe against a technician-curated Skyline result (the accepted gold standard).

THE ANSWER, IN THREE NUMBERS

R² 0.984

agreement with the curated gold standard (slope 0.998)

~8 min

to process all 918 samples vs >4 h of manual review



≤ CV

equal-or-tighter replicate precision at every level

NEXT PAGE**Agreement on standards & QCs****THEN****Calibration slope & replicate CV****FINALLY****Scientific takeaway & bottom line**

Lower limit of quantitation & cal-replicate survival

For each peptide and pipeline, the LLOQ is the lowest calibrator level whose ≥ 2 replicates survive the 25 % recovery QC. Lower is better.



LLOQ overlap

peptide	MRMpipe	Skyline	verdict
EIVPVLVSTR	0.176	0.176	= match
GPFLLGK	0.176	0.176	= match
VASHTPSLLR	0.176	0.176	= match
ITDTIGPTETSIAPR	0.176	0.176	= match
GGPLDGTYS	0.176	0.176	= match
ATISNDGATILK	0.176	0.176	= match
FLEEHPGGEEVLR	0.176	0.176	= match
AEISGVQGSETQVLR	0.176	0.176	= match
QQASEADAAAATFR	0.176	0.176	= match
ALGPEQLLR	0.176	0.176	= match
AIQDPAFSDGIR	0.176	0.176	= match
VIPSEDNLISEVEK	0.176	0.176	= match
TFITIGDR	0.176	0.176	= match
DLIEAIR	0.176	0.176	= match
ADILQVGLR	0.176	0.176	= match
LSDQFHDILIR	0.176	0.176	= match
EVVEEAENGR	0.176	0.176	= match
LSYSDSDLK	0.176	0.176	= match


peptide	MRMpipe	Skyline	verdict
VIVDFSSPNIK	0.176	0.176	= match
LNSAIYDR	0.176	0.176	= match
GADVWFK	0.176	0.176	= match
FAQNTLGK	0.176	0.176	= match
SFWSYAFSR	0.176	0.176	= match
LPSGEDYNLK	0.176	0.176	= match
VTEGLTDVILYHQDDK	0.176	0.176	= match
NELSGALTGLTR	0.176	0.176	= match
LVLDSVK	0.176	0.176	= match
IFINLPR	0.176	0.176	= match
GAEAAANTGPGGVPVQGSK	0.176	0.176	= match
GDGVLVAPPLR	0.352	0.352	= match
AGVDPLVPLR	0.352	0.352	= match
NLENGALQPSDLDR	0.352	0.352	= match
GAEAAANTGPDGVPVEGSR	0.352	0.352	= match
TLHEVLEK	7.08	7.08	= match
LFTVYPSTK	14.2	14.2	= match

Both pipelines agree on standards + QCs


1,076 standards (black) + 540 QCs (orange), above LLOQ in both pipelines.

Bottom: per-QC-level % difference (standards reserved for calibration).

Identity scatter





Per-QC-level percent difference · QC samples only





XIC — example picks

Skyline pick (left) vs MRMpipe pick (right) on the same chromatogram.

GAEAANVTGPGGVPVQGSK · STANDARD @ 202.376 ng/mL · L/H ratio SK 3.734 MP 3.729



AEISGVQGSSETQVLR · QC @ 13.136 ng/mL · L/H ratio SK 0.314 MP 0.262




Light = black, heavy IS = red. Shaded band = integration window.

Calibration slope agreement + QC replicate CV


Top: calibration slope MRMpipe vs Skyline across 35 peptides validated in both.

Bottom: replicate CV per QC concentration, all peptides pooled.

Slope vs slope · 35 peptides



QC replicate CV by concentration · all peptides pooled




Accuracy — distance from the known nominal concentration

Top: recovery box per QC level per pipeline. 100 % = perfect.


Bottom: patient unknowns (no nominal truth) — MRMpipe vs Skyline scatter.

QC recovery — distance from known truth, by pipeline



accuracy metric	MRMpipe	Skyline
median recovery	98.7 %	97.6 %
bias from 100 %	1.3	2.4
within ± 25 %	75.3 %	68.6 %
P90 $ \Delta$ from 100	37.5	41.7

Patient unknowns (n=16,779)



QC accuracy — pooled across all QC concentrations. Each pipeline's back-calculated QC vs known nominal value. High-conc is prone to fluctuations (heteroskedasticity, $1/x^2$ cal pulled by low-conc points); plus a baseline of pipetting/instrumental error common to both pipelines.

Summary — scientific takeaway

Spotlight peptide: FLEEHPGGEEVLR.

Speed.

MRMpipe processes 918 samples end-to-end in ~8 min wall time (fully automated).


Reproducibility.

MRMpipe replicate CV on QCs is tighter on 79/135 (59 %) of (peptide × level) pairs. Median CV 13.5 % vs 14.7 %.


Accuracy.

MRMpipe $|\text{median QC recovery} - 100\%| = 1.3\%$, Skyline = 2.4 %. MRMpipe slightly closer to nominal.

FLEEHPGGEEVLR · calibration curves



FLEEHPGGEEVLR · QC recovery



metric	MRMpipe	Skyline
analysis time (918 samples)	~8 minutes	manual revision
median QC recovery	98.7 %	97.6 %
median QC replicate CV	13.5 %	14.7 %

Bottom line

MRMpipe reproduces Skyline with slightly tighter CVs and better QC recovery, while being fully automated.